

med 5

NEUTRAL DISINFECTANT

MED 5 Neutral Disinfectant is a 5 minutes acting disinfectant that cleans, deodorizes and cuts through tough grease and grime. This product contains no phosphates and delivers non-acid disinfectant and cleaning performance in an economical concentrate.

It meets AOAC efficacy standards for hospital disinfectants in 5 minutes! It is effective against a broad spectrum of bacteria & viruses.

MED 5 Neutral Disinfectant has been formulated to aid in the reduction of cross-contamination on treated surfaces not only in hospitals, but in homes, schools, institutions and industry.

USES:

- Appliances
- Bathrooms
- Desks/Chairs
- Floors
- Glazed Ceramic & Glazed Porcelain
- Upholstery, vinyl & plastic
- Offices
- Sinks
- Stainless Steel
- Stoves
- Laminated Surfaces, etc.

NON-CRITICAL MEDICAL DEVICE SURFACES:

- Beds/Bed frames
- Crutches
- Examining tables
- Stethoscopes
- Resuscitators
- Rescue tools
- Defibrillators
- Medical Lamps
- Gurneys
- Stretchers
- Walkers, etc.

BENEFITS

- Neutral pH, non-staining, non-corrosive to equipment
- Easy to handle
- Does not dull floor finishes
- Dwell time: 5 minutes
- Cleaner – Bactericide – Virucide – Fungicide – Germicide – Controls Mould & Mildew
- Approved by Health Canada
- DIN 02505622

CODE & AVAILABLE FORMAT

- 3.78L INO-MED5-378



med 5

NEUTRAL DISINFECTANT

MED 5 Neutral Disinfectant is a 5 minutes acting disinfectant that cleans, deodorizes and cuts through tough grease and grime. It is effective against a broad spectrum of bacteria & viruses.



CHARACTERISTICS

Appearance	Liquid
Colour	Green
DIN	02470802
Relative Density	0.9983
Odour/perfume/fragrance	Citrus
PH	7-8
Solubility	Soluble
Viscosity at 25°C	2.024 mm ² /s@ 25° C

ACTIVE INGREDIENTS

Diocetyl Dimethyl Ammonium Chloride.....5.26%
 Alkyl (C12, 67%; C14, 25%; C16, 7%; C18, 1%)
 dimethyl benzyl ammonium chloride.....3.51%

DILUTIONS

■ 1:16 ■ 1:64

NOTE: Improper use of dilutions may cause this product not to work properly and may increase the risk of health effects

MED 5 IS EFFECTIVE AGAINST:

BACTERIA – 5 MINUTES CONTACT TIME (1:64 - 16 ml /L)

• Pseudomonas aeruginosa	• Escherichia coli – Extended Spectrum Beta Lactamase producing
• Salmonella enterica	• Escherichia coli NDM-1 - Carbapenem Resistant
• Staphylococcus aureus	• Klebsiella pneumoniae
• Acinetobacter baumannii	• Klebsiella pneumoniae – Extended Spectrum Beta Lactamase producing
• Enterobacter aerogenes	• Klebsiella pneumoniae NDM 1- Carbapenem Resistant
• Enterobacter cloacae NDM 1 • – Carbapenem Resistant	• Staphylococcus aureus - Community Associated Methicillin-Resistant (USA300)
• Enterococcus faecalis • – Vancomycin resistant	• Staphylococcus aureus - Methicillin-Resistant
• Escherichia coli	• Staphylococcus aureus - Vancomycin Resistant
• Escherichia coli O157:H7	• Streptococcus pyogenes

med 5

NEUTRAL DISINFECTANT

MED 5 Neutral Disinfectant is a 5 minutes acting disinfectant that cleans, deodorizes and cuts through tough grease and grime. It is effective against a broad spectrum of bacteria & viruses.

MED 5 IS EFFECTIVE AGAINST:

VIRUSES – 5 MINUTES CONTACT TIME (1:64 - 16 ml /L)

• Hepatitis B Virus	• Human Coronavirus
• Hepatitis C Virus	• Influenza A Virus (H7N9)
• Herpes Simplex Virus Type 1	• Influenza A virus / Hong Kong
• Herpes Simplex Virus Type 2	• Norwalk Virus – Norovirus
• HIV-1 (AIDS virus)	• Respiratory Syncytial Virus
• SARS associated Coronavirus	• Vaccinia

FUNGI

FUNGI	CONTACT TIME
• Candida albicans (1:64 -16 mL/L)	5 minutes
• Trichophyton mentagrophytes (1:16 - 64 mL/L)	10 minutes

